

E-Scope® II

Electronic Stethoscope

Improve the way you work

By improving the way you listen



Introducing a new type of stethoscope

- 🔊 **Sounds** - The maximum volume of the E-Scope II is approximately 30 times louder than an acoustic stethoscope!
- 🔊 **New Specialist Adult Diaphragm** - new stainless steel chest piece included with every E-Scope
- 🔊 **Lifetime Technical Support** - call us and a live technician answers the phone
- 🔊 **2 minute timer** - for battery conservation
- 🔊 **Models for Clinical use, EMS, hearing impaired** - the most complete product line



*E-Scope Hearing Impaired
Model for ITE hearing aids*

- 🔊 **Listen** - to the sounds through headphones
- 🔊 **Belt Clip** - swivels for ease of use
- 🔊 **Use over hearing aids** - or for comfort

Cardionics
leaders in auscultation





The E-Scope Electronic Stethoscope utilizes state of the art electronic technology to amplify heart, breath and Korotkoff sounds so they become much **easier to hear**. The sound is pure and clean. Compared to acoustic stethoscopes, murmurs and faint blows are much more distinguishable (see the below chart).

The E-Scope is **easy to use**. All controls are on the chest piece of the stethoscope. The controls include:

-  Volume Up
-  Volume Down
-  Heart or Breath Sounds switch

Afraid you might leave the E-Scope ON and lose **battery life**? The E-Scope has a built in timer that automatically shuts off the unit after 2 minutes. To reactivate the E-Scope simply press any control button on the chest piece and you're back in business. The E-Scope remembers the last volume setting upon restart. Battery life is estimated at approximately 6 months using an off-the-shelf AAA battery!

Additional features include:

-  External output for recording or telemedicine
-  3 types of eartips
-  Softpack carrying case
-  Operator's manual & instructional CD-ROM

Hearing Impaired? The E-Scope Hearing Impaired version is made just for you. Simply place the supplied pair of headphones over your in-the-ear hearing aids, or via direct audio input cable or silhouette cord for behind-the-ear hearing aids and hear heart/breath sounds with unparalleled clarity.

EMS - paramedic? The E-Scope EMS provides exceptional sound quality in a very high noise area.

For more information, uses and updates visit us at www.cardionics.com, or call us at 800-364-5901.



E-Scope EMS Electronic Stethoscope



The 'Specialist' Chest Pieces - Pediatric Bell , Adult Bell, Adult Diaphragm

The 'Specialist' heads consist of a stainless steel bell and diaphragm with a Delrin non-chill ring. They provide a 50% increase in sound sensitivity, a reduction in outside noise and comfortable feel.

Ordering Information

Cat. No.

718-7700 - Clinical Model E-Scope with Specialist adult diaphragm, chest piece/earpiece kit and carrying case, 1 year warranty

718-7740 - Same as cat no. 718-7700 except supplied with 'Specialist' Bell and Specialist Adult diaphragm

718-7710 - **Hearing Impaired** E-Scope (for ITE) with Specialist Adult Diaphragm, chest piece/ear piece kit, carrying case, headphone, 1 year warranty

718-7750 - Same as cat no. 718-7710 except supplied with 'Specialist' Adult Bell & Specialist Adult Diaphragm

718-7800 - E-Scope EMS for use in a noisy ambulance

Additional Equipment


360-7750 - Specialist Adult Bell

360-7755 - Specialist Pediatric Bell

718-0405 - External Headphones for second listener or hearing impaired

718-7030 - Pocket Monitor software for PDA or PC - displays, records & measures duration and frequency of murmurs.

E-Scope vs Acoustic Stethoscope



Cardionics, Inc.

910 Bay Star Blvd.
Telephone : 281-488-5901

Webster, Texas 77598 USA
US and Canada : 1-800-364-5901